


# Ball Head Technology®

Pioneers in Angulation



CATALOGUE  
CAD/CAM PROSTHESIS




**30° OF FREEDOM**



The evolution towards customization in dental prosthetic design increasingly requires the use of angled abutments to provide versatility, simplicity and ease of use. Our Premium and Universal Solutions are pioneers in covering these needs.



Our main goal is achieving excellence to generate safety, reliability and comfort when using our products, both for the dental lab and clinical professionals, responsible for the patient's treatment.



Since our first patent in 2008 we have developed the most efficient and innovative angled solutions based on the vision and experience of the most qualified restoration and oral surgery professionals.

# INDEX

---

|                                        |    |
|----------------------------------------|----|
| <b>1. CAD/CAM PROSTHESIS</b>           |    |
| 1.1 CAD/CAM IMPLANT LEVEL .....        | 5  |
| 1.2 CAD/CAM WELDING TO CrCo BASE ..... | 6  |
| 1.3 CAD/CAM CEMENTING TO Ti BASE ..... | 7  |
| <b>2. PRODUCTS</b> .....               | 8  |
| <b>3. BHS30 SCREWDRIVER</b> .....      | 35 |
| <b>4. PROTOCOLS</b> .....              | 36 |
| <b>5. COMPATIBILITY CHART</b> .....    | 42 |

---

Av. de Les Garrigues 17  
25001 Lleida · Spain  
**(+34) 973 831 391**  
info@ballheadtechnology.com  
www.ballheadtechnology.com

All names and trademarks used in this catalogue have been registered and are propriety of their respective companies.  
Images published are merely intended for commercial purposes and have no contractual nature.  
The information contained in this catalogue is aimed exclusively to professionals in the dental sector.  
Ball Head Technology® reserves the right to make changes to this document without prior notice.  
The guarantee conditions of Ball Head Technology® will apply as long as the instructions for each reference are respected.



## 1.1 CAD/CAM IMPLANT LEVEL

---



### ABSOLUTE FREEDOM AND ACCESSIBILITY

---

Wide flat wall for an optimum Best Fit (Position and Orientation).

Manufactured in Titanium Grade 5 with tolerances below  $\pm 5\mu\text{m}$ .

Optimal surface finishing for a perfect scanning experience.

Exact connections for all major implant systems.

## 1.2 CAD/CAM WELDING TO CrCo BASE

---



### AFFORDABLE STRUCTURES WITH PERFECT MECHANIZED SEATING

---

#### **CrCo Base:**

Optimum height allowing 30° angulation for the BHS screwdriver.

Conical design for multiple restorations and engaging flat wall for single ones.

Provides an easy passive fit in CrCo sinterized multiple restorations.

#### **Scanbody:**

Wide flat wall for an optimum Best Fit (Position and Orientation).

Manufactured in Titanium Grade 5 with tolerances below  $\pm 5\mu\text{m}$ .

Optimal surface finishing for a perfect scanning experience.

Exact connections for all major implant systems.

### 1.3 CAD/CAM CEMENTING TO Ti BASE

---



#### PROVIDING AESTHETIC AND FUNCTIONAL RESULTS

---

##### **Ti Base:**

U-shape cut for angled screwdriver access up to 30°.

Double vestibular height for a safe cementation.

Indented body design to improve axial retention.

Optimum emergence profile for soft tissue care. Ask for gingival heights

##### **Scanbody:**




Wide flat wall for an optimum Best Fit (Position and Orientation).

Manufactured in Titanium Grade 5 with tolerances below  $\pm 5\mu\text{m}$ .


Optimal surface finishing for a perfect scanning experience.

Exact connections for all major implant systems.



## 2. PRODUCTS

|                                                                                     | PLATFORM Ø          | SCANBODY | LIBRARY PLATFORM    |                                                                                           | CrCo BASE            | Ti BASE              | CLINICAL SCREW | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|---------------------|----------|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|----------------------|----------------|-----------|-----------------|
| AB IMPLANTS COMPATIBLE                                                              |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
| I SERIES INTERNAL HEX                                                               |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
|    | Narrow Ø 3 - Ø 3,75 | BHS08951 | BHS30-9500 (NARROW) | <br>     | BHS03200<br>BHS03201 | BHS03198<br>BHS03199 | BHS02041       | BHS02042  | EM045130        |
|                                                                                     | STD Ø 3,5 - 6,0     | BHS08321 | BHS30-3200 (NP)     | <br>     | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  |                 |
| I 10C CONUS                                                                         |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
|    | RP Ø 4,2            | BHS08312 | BHS30 - 3100 (RP)   | <br>     | BHS03046<br>BHS03047 | BHS03044<br>BHS03045 | BHS02023       | BHS02024  | EM048149        |
| ACE COMPATIBLE                                                                      |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
| INFINITY EXTERNAL HEX                                                               |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
|    | NP Ø 3,3 / 3,75     | BHS08611 | BHS30 - 6100 (NP)   | <br>     |                      | BHS03010<br>BHS03011 | BHS02005       | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,0            | BHS08612 | BHS30 - 6100 (RP)   | <br>     | BHS03003<br>BHS03004 | BHS03005<br>BHS03006 | BHS02001       | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0            | BHS08613 | BHS30 - 6100 (WP)   |                                                                                                                                                                            |                      |                      | BHS02001       | BHS02002  | EM045082        |
| INFINITY TRI - CAM                                                                  |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | NP Ø 3,5            | BHS08111 | BHS30 - 1100 (NP)   | <br> |                      | BHS03060<br>BHS03061 | BHS02027       | BHS02028  | EM046291        |
|                                                                                     | RP Ø 4,3            | BHS08112 | BHS30 - 1100 (RP)   | <br> |                      | BHS03012<br>BHS03013 | BHS02009       | BHS02010  | EM046292        |
|                                                                                     | WP Ø 5,0            | BHS08113 | BHS30 - 1100 (WP)   |                                                                                                                                                                            |                      |                      | BHS02009       | BHS02010  | EM046293        |
| INFINITY INTERNAL HEX                                                               |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | NP Ø 3,7 / 4,1      | BHS08321 | BHS30 - 3200 (NP)   | <br> | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  | EM045130        |
|                                                                                     | RP Ø 4,7 / 5,1      | BHS08322 | BHS30 - 3200 (RP)   | <br> | BHS03033<br>BHS03034 | BHS03020<br>BHS03021 | BHS02021       | BHS02022  | EM046218        |
| ADIN COMPATIBLE                                                                     |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
| SWELL                                                                               |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | Ø 3,5 - 6           | BHS08321 | BHS30 - 3200 (NP)   | <br> | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  | EM045130        |
| TOUAREG S / OS                                                                      |                     |          |                     |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | Ø 3,5 - 6           | BHS08321 | BHS30 - 3200 (NP)   | <br> | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  | EM045130        |


## 2. PRODUCTS

|                                                                                     | PLATFORM Ø               | SCANBODY | LIBRARY PLATFORM         |                                                                                           | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|--------------------------|----------|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>TOUAREG CLOSEFIT</b>                                                             |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RP Ø 3,5                 | BHS08311 | <b>BHS30 - 3100 (NP)</b> | <br>     | <b>BHS03058</b><br><b>BHS03059</b> | <b>BHS03056</b><br><b>BHS03057</b> | <b>BHS02025</b> | BHS02026  | EM048148        |
|                                                                                     | WP Ø 4,3 / 5             | BHS08312 | <b>BHS30 - 3100 (RP)</b> |                                                                                                                                                                            | <b>BHS03046</b><br><b>BHS03047</b> | <b>BHS03044</b><br><b>BHS03045</b> | <b>BHS02023</b> | BHS02024  | EM048149        |
| <b>ALPHABIO COMPATIBLE</b>                                                          |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EXTERNAL ALPHABIO</b>                                                            |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5                 | BHS08211 | <b>BHS30 - 2100 (NP)</b> | <br>     |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1                 | BHS08212 | <b>BHS30 - 2100 (RP)</b> | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1                 | BHS08213 | <b>BHS30 - 2100 (WP)</b> |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>INTERNAL ALPHABIO</b>                                                            |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,3 - 6                | BHS08321 | <b>BHS30 - 3200 (NP)</b> | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>ALPHA DENT COMPATIBLE</b>                                                        |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>ACTIVE</b>                                                                       |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,75 / 4,2 / 5,0       | BHS08321 | <b>BHS30 - 3200 (NP)</b> | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>ACTIVE PLUS</b>                                                                  |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 4,2 / 5,0              | BHS08321 | <b>BHS30 - 3200 (NP)</b> | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>CLASSIC</b>                                                                      |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,3 / 3,75 / 4,2 / 5,0 | BHS08321 | <b>BHS30 - 3200 (NP)</b> | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>ANCLADEN COMPATIBLE</b>                                                          |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>ANCLALOCK</b>                                                                    |                          |          |                          |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,75 - 5               | BHS08321 | <b>BHS30 - 3200 (NP)</b> | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |


## 2. PRODUCTS

|                                                                                     | PLATFORM Ø   | SCANBODY | LIBRARY PLATFORM                       | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|--------------|----------|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>ANTHOGRYR COMPATIBLE</b>                                                         |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>ANTHOFIT EH</b>                                                                  |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RP Ø 4,1     | BHS08612 | <b>BHS30 - 6100 (RP)</b>               | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0     | BHS08613 | <b>BHS30 - 6100 (WP)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>OSSFIT</b>                                                                       |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RP Ø 4,8     | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b>         | <br>     |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
|                                                                                     | SYN_RP Ø 4,8 | BHS08713 | <b>BHS30 - 7100 (RN SYNOCTA Ø 4,8)</b> |                                                                                           |                                    | <b>BHS03019</b>                    | <b>BHS02013</b> | BHS02014  | EM046375        |
| <b>ARDS COMPATIBLE</b>                                                              |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>SMART</b>                                                                        |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,75 - 4,5 | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>CLASSIC</b>                                                                      |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,3 - 6    | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>PREMIUM</b>                                                                      |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,3 - 6    | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>ASTRA COMPATIBLE</b>                                                             |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>AQUA</b>                                                                         |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RP_3.5 / 4.0 | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 / 4.0)</b>   | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
| <b>LILAC</b>                                                                        |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | WP_4.5 / 5.0 | BHS08123 | <b>BHS30 - 1200 (WP Ø 4.5 / 5.0)</b>   |                                                                                                                                                                            |                                    |                                    | <b>BHS02009</b> | BHS02010  | EM046370        |
| <b>CONE 20°</b>                                                                     |              |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,5        | BHS08124 | <b>BHS30 - 1200 (20°)</b>              |                                                                                         |                                    | <b>BHS03072</b>                    | <b>BHS02003</b> | BHS02004  | EM045638        |

## 2. PRODUCTOS

|                                                                                     | PLATFORM Ø                  | SCANBODY | LIBRARY PLATFORM           | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|-----------------------------|----------|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>CONE 45°</b>                                                                     |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,5                       | BHS08125 | <b>BHS30 - 1200 (45°)</b>  |                                                                                           |                                    | <b>BHS03074</b>                    | <b>BHS02003</b> | BHS02004  | EM047786        |
|                                                                                     |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EV</b>                                                                           |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø3,6                        | BHS08152 | <b>BHS30 - 1500 (Ø3,6)</b> | <br>     | <b>BHS03154</b><br><b>BHS03155</b> | <b>BHS03152</b><br><b>BHS03153</b> | <b>BHS02053</b> | BHS02054  |                 |
|                                                                                     | Ø4,2                        | BHS08153 | <b>BHS30 - 1500 (Ø4,2)</b> | <br>     | <b>BHS03158</b><br><b>BHS03159</b> | <b>BHS03156</b><br><b>BHS03157</b> | <b>BHS02055</b> | BHS02056  |                 |
| <b>AVINENT COMPATIBLE</b>                                                           |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EH</b>                                                                           |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5                    | BHS08211 | <b>BHS30 - 2100 (NP)</b>   | <br>     |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02031</b> | BHS02032  | EM045222        |
|                                                                                     | RP Ø 4,1                    | BHS08212 | <b>BHS30 - 2100 (RP)</b>   | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1                    | BHS08213 | <b>BHS30 - 2100 (WP)</b>   |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>IH</b>                                                                           |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,5 / 4,1                 | BHS08321 | <b>BHS30 - 3200 (NP)</b>   | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>MULTIUNIT</b>                                                                    |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RP Ø 4,8                    | BHS08411 | <b>BHS30 - 4100 (RP)</b>   |                                                                                         | <b>BHS03008</b>                    | <b>BHS03009</b>                    | <b>BHS02003</b> | BHS02004  | EM045169        |
| <b>OCEAN CONICAL (C.C.)</b>                                                         |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,5 (Implant Ø 3,5/Ø 4,0) | BHS08311 | <b>BHS30 - 3100 (NP)</b>   | <br> | <b>BHS03058</b><br><b>BHS03059</b> | <b>BHS03056</b><br><b>BHS03057</b> | <b>BHS02025</b> | BHS02026  | EM048148        |
|                                                                                     | Ø 4,1 (Implant Ø 4,5/Ø 5,0) | BHS08312 | <b>BHS30 - 3100 (RP)</b>   | <br> | <b>BHS03046</b><br><b>BHS03047</b> | <b>BHS03044</b><br><b>BHS03045</b> | <b>BHS02023</b> | BHS02024  | EM048149        |
| <b>BEGO COMPATIBLE</b>                                                              |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EH</b>                                                                           |                             |          |                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP_Ø 2,7 / 3,1              | BHS08421 | <b>BHS30 - 4200 (NP)</b>   |                                                                                                                                                                            |                                    |                                    | <b>BHS02013</b> | BHS02014  | EM047957        |



## 2. PRODUCTS

|                                                                                     | PLATFORM Ø           | SCANBODY | LIBRARY PLATFORM                       | <br>     | CrCo BASE                          | Ti BASE                                                                           | CLINICAL SCREW  | LAB SCREW       | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|----------------------|----------|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|-----------------------------------------------------------------------------------|-----------------|-----------------|-----------------|
| IH                                                                                  |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
|    | Ø 3,75 / 4,1         | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b>                                                | <b>BHS02021</b> | BHS02022        | EM045130        |
| <b>BIOGENESIS COMPATIBLE</b>                                                        |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
| 3ICON                                                                               |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
|    | NP Ø 3,5             | BHS08611 | <b>BHS30 - 6100 (NP)</b>               | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b>                                                | <b>BHS02001</b> | BHS02002        | EM045468        |
|                                                                                     | RP Ø 4,1             | BHS08612 | <b>BHS30 - 6100 (RP)</b>               |                                                                                                                                                                            |                                    |                                                                                   | <b>BHS02001</b> | BHS02002        | EM045098        |
|                                                                                     | WP Ø 5,0             | BHS08613 | <b>BHS30 - 6100 (WP)</b>               |                                                                                                                                                                            |                                    |                                                                                   | <b>BHS02001</b> | BHS02002        | EM045082        |
| ITICON                                                                              |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
|    | RP Ø 4,8             | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b>         | <br>     |                                    | <b>BHS03016</b><br><b>BHS03017</b>                                                | <b>BHS02011</b> | BHS02012        | EM045144        |
|    | RP Ø 4,8             | BHS08713 | <b>BHS30 - 7100 (RN SYNOCTA Ø 4,8)</b> |                                                                                                                                                                            |                                    |  | <b>BHS03019</b> | <b>BHS02013</b> | BHS02014        |
| ATICON                                                                              |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
|  | RP Ø 3,5 / 4,0       | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 / 4.0)</b>   | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b>                                                | <b>BHS02025</b> | BHS02026        | EM046371        |
|                                                                                     | WP Ø 4,5 / 5,0       | BHS08123 | <b>BHS30 - 1200 (WP Ø 4.5 / 5.0)</b>   |                                                                                                                                                                            |                                    |                                                                                   | <b>BHS02009</b> | BHS02010        | EM046370        |
| KISSES                                                                              |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
|  | RP Ø 3,5 / 4,0       | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 / 4.0)</b>   | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b>                                                | <b>BHS02025</b> | BHS02026        | EM046371        |
|                                                                                     | WP Ø 4,5 / 5,0 / 5,5 | BHS08123 | <b>BHS30 - 1200 (WP Ø 4.5 / 5.0)</b>   |                                                                                                                                                                            |                                    |                                                                                   | <b>BHS02009</b> | BHS02010        | EM046370        |
| CONE 20°                                                                            |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
|  | RP / WP              | BHS08124 | <b>BHS30 - 1200 (20°)</b>              |                                                                                         |                                    | <b>BHS03072</b><br><b>BHS02003</b>                                                |                 | BHS02004        | EM045638        |
| CONE 45°                                                                            |                      |          |                                        |                                                                                                                                                                            |                                    |                                                                                   |                 |                 |                 |
|  | RP / WP              | BHS08125 | <b>BHS30 - 1200 (45°)</b>              |                                                                                         |                                    | <b>BHS03074</b><br><b>BHS02003</b>                                                |                 | BHS02004        | EM047786        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø | SCANBODY | LIBRARY PLATFORM         |    | CrCo BASE       | Ti BASE         | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|------------|----------|--------------------------|-------------------------------------------------------------------------------------|-----------------|-----------------|-----------------|-----------|-----------------|
| <b>BIOHORIZONS COMPATIBLE</b>                                                       |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
| <b>INTERNAL</b>                                                                     |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
|    | NP Ø 3,5   | BHS08321 | <b>BHS30 - 3200 (NP)</b> |    | <b>BHS03050</b> | <b>BHS03048</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | RP Ø 4,5   | BHS08322 | <b>BHS30 - 3200 (RP)</b> |    | <b>BHS03051</b> | <b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
|                                                                                     | WP Ø 5,7   | BHS08323 | <b>BHS30 - 3200 (WP)</b> |    | <b>BHS03033</b> | <b>BHS03020</b> | <b>BHS02021</b> | BHS02022  |                 |
|                                                                                     |            |          |                          |    | <b>BHS03034</b> | <b>BHS03021</b> |                 |           |                 |
| <b>BIOMET 3i COMPATIBLE</b>                                                         |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
| <b>EH OSSEOTITE</b>                                                                 |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
|    | NP Ø 3,5   | BHS08611 | <b>BHS30 - 6100 (NP)</b> |                                                                                     |                 |                 | <b>BHS02001</b> | BHS02002  | EM045468        |
|                                                                                     | RP Ø 4,1   | BHS08612 | <b>BHS30 - 6100 (RP)</b> |    | <b>BHS03003</b> | <b>BHS03005</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0   | BHS08613 | <b>BHS30 - 6100 (WP)</b> |    | <b>BHS03004</b> | <b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>CERTAIN</b>                                                                      |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
|  | NP Ø 3,4   | BHS08511 | <b>BHS30 - 5100 (NP)</b> |                                                                                     |                 |                 | <b>BHS02007</b> | BHS02008  | EM045166        |
|                                                                                     | RP Ø 4,1   | BHS08512 | <b>BHS30 - 5100 (RP)</b> |                                                                                     |                 |                 | <b>BHS02007</b> | BHS02008  | EM045165        |
|                                                                                     | WP Ø 5,0   | BHS08513 | <b>BHS30 - 5100 (WP)</b> |                                                                                     |                 |                 | <b>BHS02007</b> | BHS02008  | EM045161        |
| <b>BIONER COMPATIBLE</b>                                                            |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
| <b>IKELT / BIKELT</b>                                                               |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
|  | RP Ø 4,1   | BHS08212 | <b>BHS30 - 2100 (RP)</b> |  | <b>BHS03003</b> | <b>BHS03005</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0   | BHS08613 | <b>BHS30 - 6100 (WP)</b> |  | <b>BHS03004</b> | <b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045082        |
|                                                                                     | WP Ø 5,1   | BHS08213 | <b>BHS30 - 2100 (WP)</b> |                                                                                     |                 |                 | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>HIKELT</b>                                                                       |            |          |                          |                                                                                     |                 |                 |                 |           |                 |
|  | NP Ø 3,8   | BHS08321 | <b>BHS30 - 3200 (NP)</b> |  | <b>BHS03050</b> | <b>BHS03048</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | RP Ø 4,7   | BHS08322 | <b>BHS30 - 3200 (RP)</b> |  | <b>BHS03051</b> | <b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
|                                                                                     |            |          |                          |  | <b>BHS03033</b> | <b>BHS03020</b> | <b>BHS02021</b> |           |                 |
|                                                                                     |            |          |                          |  | <b>BHS03034</b> | <b>BHS03021</b> |                 |           |                 |


## 2. PRODUCTS

|                                                                                     |                | PLATFORM Ø                    | SCANBODY                 | LIBRARY PLATFORM                                                                                                                                                           | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|----------------|-------------------------------|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| TOP DM                                                                              |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    |                | Ø 3,5 / 4,0 / 5,0             | BHS08311                 | <b>BHS30 - 3100 (NP)</b>                                                                                                                                                   | <br>     | <b>BHS03058</b><br><b>BHS03059</b> | <b>BHS03056</b><br><b>BHS03057</b> | <b>BHS02025</b> | BHS02026  | EM048148        |
| MULTIUNIT                                                                           |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    |                | RP Ø 4,8                      | BHS08411                 | <b>BHS30 - 4100 (RP)</b>                                                                                                                                                   |                                                                                           | <b>BHS03008</b>                    | <b>BHS03009</b>                    | <b>BHS02003</b> | BHS02004  | EM045169        |
| BTI COMPATIBLE                                                                      |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| EH                                                                                  |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Tiny Ø 3,5  | BHS08421                      | <b>BHS30 - 4200 (NP)</b> |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    | <b>BHS02013</b> | BHS02014  | EM047957        |
|                                                                                     | RP Ø 4,1       | BHS08422                      | <b>BHS30 - 4200 (RP)</b> | <br>     | <b>BHS03003</b><br><b>BHS03004</b>                                                                                                                                         | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b>                    | BHS02002        | EM045098  |                 |
|                                                                                     | WP Ø 5,5       | BHS08423                      | <b>BHS30 - 4200 (WP)</b> |                                                                                                                                                                            |                                                                                                                                                                            |                                    | <b>BHS02013</b>                    | BHS02014        | EM047967  |                 |
| CUATRILOBULAR                                                                       |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | RP Ø 4,1       | BHS08231                      | <b>BHS30 - 2300 (RP)</b> |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    | <b>BHS02021</b> | BHS02022  | EM047743        |
|                                                                                     | WP Ø 5,5       | BHS08232                      | <b>BHS30 - 2300 (WP)</b> |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    | <b>BHS02021</b> | BHS02022  |                 |
| BTK COMPATIBLE                                                                      |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| 8NECK INT                                                                           |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  |                | SR Ø 4,8                      | BHS08712                 | <b>BHS30 - 7100 (RN Ø 4,8 )</b>                                                                                                                                            | <br> |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
| EVO DL EXT                                                                          |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP Ø 3,4 / 3,7 | BHS08611                      | <b>BHS30 - 6100 (NP)</b> |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    | <b>BHS02035</b> | BHS02036  |                 |
|                                                                                     | RP Ø 4,1 / 4,6 | BHS08612                      | <b>BHS30 - 6100 (RP)</b> | <br> | <b>BHS03003</b><br><b>BHS03004</b>                                                                                                                                         | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b>                    | BHS02002        | EM045098  |                 |
|                                                                                     | WP Ø 5,0 / 5,6 | BHS08613                      | <b>BHS30 - 6100 (WP)</b> |                                                                                                                                                                            |                                                                                                                                                                            |                                    | <b>BHS02001</b>                    | BHS02002        | EM045082  |                 |
| EVO DL INT                                                                          |                |                               |                          |                                                                                                                                                                            |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  |                | Ø 3,5 / 3,7 / 4,6 / 5,0 / 5,6 | BHS08321                 | <b>BHS30 - 3200 (NP)</b>                                                                                                                                                   | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø       | SCANBODY | LIBRARY PLATFORM               | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|------------------|----------|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>KLASSIC EXT</b>                                                                  |                  |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,4         | BHS08611 | <b>BHS30 - 6100 (NP)</b>       |                                                                                                                                                                            |                                    |                                    | <b>BHS02035</b> | BHS02036  |                 |
|                                                                                     | RP Ø 4,1         | BHS08612 | <b>BHS30 - 6100 (RP)</b>       | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0         | BHS08613 | <b>BHS30 - 6100 (WP)</b>       |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>KLASSIC INT</b>                                                                  |                  |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,5/ 4,5 / 5,0 | BHS08321 | <b>BHS30 - 3200 (NP)</b>       | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>KONIK EXT</b>                                                                    |                  |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,4         | BHS08611 | <b>BHS30 - 6100 (NP)</b>       |                                                                                                                                                                            |                                    |                                    | <b>BHS02035</b> | BHS02036  |                 |
|                                                                                     | RP Ø 4,1         | BHS08612 | <b>BHS30 - 6100 (RP)</b>       | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0         | BHS08613 | <b>BHS30 - 6100 (WP)</b>       |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>KONIC INT</b>                                                                    |                  |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | Ø 3,5/ 5,0       | BHS08321 | <b>BHS30 - 3200 (NP)</b>       | <br>  | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>SAFETL INT</b>                                                                   |                  |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | SR Ø 4,8         | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b> | <br> |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
| <b>B&amp;W COMPATIBLE</b>                                                           |                  |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>HEXAGONO EXT</b>                                                                 |                  |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP Ø 3,3         | BHS08741 | <b>BHS30 - 7400 (S2)</b>       |                                                                                                                                                                            |                                    |                                    | <b>BHS02035</b> | BHS02036  | EM048022        |
|                                                                                     | RP Ø 4,1         | BHS08742 | <b>BHS30 - 7400 (S3 y S4)</b>  | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0         | BHS08613 | <b>BHS30 - 6100 (WP)</b>       |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø               | SCANBODY | LIBRARY PLATFORM       |                                                                                           | CrCo BASE            | Ti BASE              | CLINICAL SCREW | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|--------------------------|----------|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|----------------------|----------------|-----------|-----------------|
| HEXAGONO INT                                                                        |                          |          |                        |                                                                                                                                                                            |                      |                      |                |           |                 |
|    | NP Ø 3,3                 | BHS08321 | BHS30 - 3200 (NP)      | <br>     | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  | EM045130        |
|                                                                                     | RP Ø 4,0                 | BHS08321 | BHS30 - 3200 (NP)      | <br>     | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  | EM045130        |
|                                                                                     | WP Ø 5,0                 | BHS08321 | BHS30 - 3200 (NP)      | <br>     | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  | EM045130        |
| CAMLOG COMPATIBLE                                                                   |                          |          |                        |                                                                                                                                                                            |                      |                      |                |           |                 |
| CAMLOG                                                                              |                          |          |                        |                                                                                                                                                                            |                      |                      |                |           |                 |
|   | Ø 3,3                    | BHS08721 | BHS30-7200 (Ø3,3)      | <br>     | BHS03130<br>BHS03131 | BHS03128<br>BHS03129 | BHS02049       | BHS02050  |                 |
|                                                                                     | Ø 3,8                    | BHS08722 | BHS30-7200 (Ø3,8)      | <br>     | BHS03134<br>BHS03135 | BHS03132<br>BHS03133 | BHS02049       | BHS02050  |                 |
|                                                                                     | Ø 4,3                    | BHS08723 | BHS30-7200 (Ø4,3)      | <br>     | BHS03138<br>BHS03139 | BHS03136<br>BHS03137 | BHS02049       | BHS02050  |                 |
|                                                                                     | Ø 5,0                    | BHS08724 | BHS30-7200 (Ø5,0)      | <br>   | BHS03142<br>BHS03143 | BHS03140<br>BHS03141 | BHS02051       | BHS02052  |                 |
| CONEXÃO SIST. PROTESE COMPATIBLE                                                    |                          |          |                        |                                                                                                                                                                            |                      |                      |                |           |                 |
| MORSE INDEXADO                                                                      |                          |          |                        |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | Ø 3,5 / 3,75 / 4,0 / 5,0 | BHS08311 | BHS30 - 3100 (NP)      | <br> | BHS03058<br>BHS03059 | BHS03056<br>BHS03057 | BHS02025       | BHS02026  | EM048148        |
| HEXAGONO EXT                                                                        |                          |          |                        |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | Ø 3,3                    | BHS08741 | BHS30 - 7400 (S2)      |                                                                                                                                                                            |                      |                      | BHS02035       | BHS02036  | EM048022        |
|                                                                                     | Ø 3,75 / 4,0             | BHS08742 | BHS30 - 7400 (S3 y S4) | <br> | BHS03003<br>BHS03004 | BHS03005<br>BHS03006 | BHS02001       | BHS02002  | EM045098        |
|                                                                                     | Ø 5,0                    | BHS08613 | BHS30 - 6100 (WP)      |                                                                                                                                                                            |                      |                      | BHS02001       | BHS02002  | EM045082        |
| HEXAGONO INT                                                                        |                          |          |                        |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | Ø 3,75 / 4,0 / 5,0       | BHS08322 | BHS30 - 3200 (RP)      | <br> | BHS03033<br>BHS03034 | BHS03020<br>BHS03021 | BHS02021       | BHS02022  | EM046218        |



## 2. PRODUCTS

|                                                                                                                                                                         | PLATFORM Ø              | SCANBODY | LIBRARY PLATFORM               | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|----------|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>CORTEX COMPATIBLE</b>                                                                                                                                                |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>DYNAMIX</b>                                                                                                                                                          |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|       | Ø 3,3 / 3,8 / 4,2 / 5   | BHS08321 | <b>BHS30 - 3200 (NP)</b>       | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>CLASSIX</b>                                                                                                                                                          |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|       | Ø 3,3 / 3,8 / 4,2 / 5   | BHS08321 | <b>BHS30 - 3200 (NP)</b>       | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>SATURN</b>                                                                                                                                                           |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|       | Ø 3,8 / 4,2             | BHS08321 | <b>BHS30 - 3200 (NP)</b>       | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>COWELLMEDI COMPATIBLE</b>                                                                                                                                            |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EXTERNAL</b>                                                                                                                                                         |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|     | RP Ø 4,1                | BHS08212 | <b>BHS30 - 2600 (RP)</b>       | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                                                                                                         | WP Ø 5,1                | BHS08213 | <b>BHS30 - 2600 (WP)</b>       |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>INTERNAL</b>                                                                                                                                                         |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | RN Ø 4,8                | BHS08712 | <b>BHS30 - 7600 (RN Ø 4,8)</b> | <br> |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
| <b>SUBMERGED (INT)</b>                                                                                                                                                  |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | Ø 3,5 / 4 / 4,5 / 5 / 6 | BHS08432 | <b>BHS30 - 4600 (4.0)</b>      | <br> | <b>BHS03102</b><br><b>BHS03103</b> | <b>BHS03100</b><br><b>BHS03101</b> | <b>BHS02009</b> | BHS02010  |                 |
| <b>DENTALTECH COMPATIBLE</b>                                                                                                                                            |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>IMPLOGIC GII</b>                                                                                                                                                     |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | RP Ø 4,5                | BHS08322 | <b>BHS30 - 3200 (RP)</b>       | <br> | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
| <b>IMPLASSIC HX2</b>                                                                                                                                                    |                         |          |                                |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | RP Ø 4,1                | BHS08212 | <b>BHS30 - 2100 (RP)</b>       | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |

## 2. PRODUCTS

|                                                                                                                                                                         |                               | PLATFORM Ø | SCANBODY                | LIBRARY PLATFORM                                                                    | <br> | CrCo BASE | Ti BASE  | CLINICAL SCREW | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|------------|-------------------------|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|----------|----------------|-----------|-----------------|
| DENTIS COMPATIBLE                                                                                                                                                       |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
| S-CLEAN                                                                                                                                                                 |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
|       | RP Ø 3,4                      | BHS08432   | BHS30 - 4600 (4.0)      |    | BHS03102                                                                                                                                                               | BHS03100  | BHS02009 | BHS02010       |           |                 |
|                                                                                                                                                                         |                               |            |                         |    | BHS03103                                                                                                                                                               | BHS03101  |          |                |           |                 |
|       | RP Ø 4,8                      | BHS08712   | BHS30 - 7100 (RN Ø 4,8) |    |                                                                                                                                                                        | BHS03016  | BHS02011 | BHS02012       |           | EM045144        |
|                                                                                                                                                                         |                               |            |                         |    |                                                                                                                                                                        | BHS03017  |          |                |           |                 |
| DENTIUM COMPATIBLE                                                                                                                                                      |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
| SIMPLE LINE II                                                                                                                                                          |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
|       | RN Ø 4,8                      | BHS08712   | BHS30 - 7600 (RN Ø 4,8) |    |                                                                                                                                                                        | BHS03016  | BHS02011 | BHS02012       |           | EM045144        |
|                                                                                                                                                                         |                               |            |                         |    |                                                                                                                                                                        | BHS03017  |          |                |           |                 |
| IMPLANTUM                                                                                                                                                               |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
|   | Ø 3,6 / 4 / 4,4 / 4,9 / 6 / 7 | BHS08432   | BHS30 - 4600 (4.0)      |  | BHS03102                                                                                                                                                               | BHS03100  | BHS02009 | BHS02010       |           |                 |
|                                                                                                                                                                         |                               |            |                         |  | BHS03103                                                                                                                                                               | BHS03101  |          |                |           |                 |
| SUPER LINE                                                                                                                                                              |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
|   | Ø 3,6 / 4 / 4,5 / 5 / 6       | BHS08432   | BHS30 - 4600 (4.0)      |  | BHS03102                                                                                                                                                               | BHS03100  | BHS02009 | BHS02010       |           |                 |
|                                                                                                                                                                         |                               |            |                         |  | BHS03103                                                                                                                                                               | BHS03101  |          |                |           |                 |
| DENTSPLY SIRONA COMPATIBLE                                                                                                                                              |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
| XIVE S                                                                                                                                                                  |                               |            |                         |                                                                                     |                                                                                                                                                                        |           |          |                |           |                 |
|                                                                                      | D 3.4                         | BHS08221   | BHS30-2200 (D3.4)       |  | BHS03118                                                                                                                                                               | BHS03116  | BHS02047 | BHS02048       |           |                 |
|                                                                                                                                                                         |                               |            |                         |  | BHS03119                                                                                                                                                               | BHS03117  |          |                |           |                 |
|                                                                                                                                                                         | D 3.8                         | BHS08222   | BHS30-2200 (D3.8)       |  | BHS03122                                                                                                                                                               | BHS03120  | BHS02047 | BHS02048       |           |                 |
|                                                                                                                                                                         |                               |            |                         |  | BHS03123                                                                                                                                                               | BHS03121  |          |                |           |                 |
|                                                                                                                                                                         | D 4.5                         | BHS08223   | BHS30-2200 (D4.5)       |  | BHS03126                                                                                                                                                               | BHS03124  | BHS02047 | BHS02048       |           |                 |
|                                                                                                                                                                         |                               |            |                         |  | BHS03127                                                                                                                                                               | BHS03125  |          |                |           |                 |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø                     | SCANBODY | LIBRARY PLATFORM                    |                                                                                           | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|--------------------------------|----------|-------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>DIO IMPLANTS COMPATIBLE</b>                                                      |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EXTERNAL</b>                                                                     |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5                       | BHS08211 | <b>BHS30 - 2100 (NP)</b>            |                                                                                           |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1                       | BHS08212 | <b>BHS30 - 2100 (RP)</b>            | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1                       | BHS08213 | <b>BHS30 - 2100 (WP)</b>            |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>UF II</b>                                                                        |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,8 / 4 / 4,5 / 5 / 5,5      | BHS08432 | <b>BHS30 - 4600 (4.0)</b>           | <br>     | <b>BHS03102</b><br><b>BHS03103</b> | <b>BHS03100</b><br><b>BHS03101</b> | <b>BHS02009</b> | BHS02010  |                 |
| <b>EDISON MEDICAL COMPATIBLE</b>                                                    |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>LAMINA (IH)</b>                                                                  |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,3 / 3,75 / 4,2 / 5,0 / 6,0 | BHS08321 | <b>BHS30 - 3200 (NP)</b>            | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>FORTIS (IH)</b>                                                                  |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | Ø 3,3 / 3,75 / 4,2 / 5,0 / 6,0 | BHS08321 | <b>BHS30 - 3200 (NP)</b>            | <br>   | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>ELITE MEDICA COMPATIBLE</b>                                                      |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EXTERNAL CONNECT</b>                                                             |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP Ø 3,5                       | BHS08211 | <b>BHS30 - 2100 (NP)</b>            |                                                                                         |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1                       | BHS08212 | <b>BHS30 - 2100 (RP)</b>            | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1                       | BHS08213 | <b>BHS30 - 2100 (WP)</b>            |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>EUROTEKNIKA COMPATIBLE</b>                                                       |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>UNEVA</b>                                                                        |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,6 / 4,1 / 4,8              | BHS08212 | <b>BH30 - 2100 (RP)</b>             | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
| <b>NATURACTIS</b>                                                                   |                                |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,5 / 4 / 4,5 / 5            | BHS08122 | <b>BH30 - 1200 (RP Ø 3.5 / 4.0)</b> | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø             | SCANBODY | LIBRARY PLATFORM                        | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|------------------------|----------|-----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>NATURALL +</b>                                                                   |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,5 / 4 / 4,5 / 5    | BHS08122 | <b>BHS30 - 1200</b><br>(RP Ø 3.5 / 4.0) | <br>     | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
|                                                                                     |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>NATEA</b>                                                                        |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,6 / 4,1         | BHS08122 | <b>BHS30 - 1200</b><br>(RP Ø 3.5 / 4.0) | <br>     | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
|                                                                                     | RP Ø 4,8 / 6           | BHS08122 | <b>BHS30 - 1200</b><br>(RP Ø 3.5 / 4.0) | <br>     | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
| <b>AESTHETICA</b>                                                                   |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Ø 3,6 / 4,1 / 4,8      | BHS08712 | <b>BHS30 - 7100</b><br>(RN Ø 4,8)       | <br>     |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
| <b>GALIMPLANT COMPATIBLE</b>                                                        |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>IPXE - EXTERNA</b>                                                               |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RP Ø 3,5 / 4 / 5       | BHS08212 | <b>BHS30 - 2100</b><br>(RP)             | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
| <b>ICI - INTERNA</b>                                                                |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RP Ø 3,2 / 3,5 / 4 / 5 | BHS08122 | <b>BHS30 - 1200</b><br>(RP Ø 3.5 / 4.0) | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
| <b>IPX - INTERNA</b>                                                                |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RP Ø 3,5 / 4 / 4,5 / 5 | BHS08122 | <b>BHS30 - 1200</b><br>(RP Ø 3.5 / 4.0) | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
| <b>GMI COMPATIBLE</b>                                                               |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>FRONTIER (IH)</b>                                                                |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RP Standard Ø 3,3      | BHS08321 | <b>BHS30 - 3200</b><br>(NP)             | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | WP Wide Ø 4,3          | BHS08322 | <b>BHS30 - 3200</b><br>(RP)             | <br> | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
| <b>FRONTIER multi-aesthetic</b>                                                     |                        |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RP Ø 4,8               | BHS08411 | <b>BHS30 - 4100</b><br>(RP)             |                                                                                         | <b>BHS03008</b>                    | <b>BHS03009</b>                    | <b>BHS02003</b> | BHS02004  | EM045169        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø            | SCANBODY | LIBRARY PLATFORM                       | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|-----------------------|----------|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>PHOENIX (EH)</b>                                                                 |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Standard Ø 4,1        | BHS08212 | <b>BHS30 - 2100 (RP)</b>               | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | Wide (Hex. 2,7) Ø 5,1 | BHS08613 | <b>BHS30 - 6100 (WP)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
|                                                                                     | Wide (Hex. 3,4) Ø 5,1 | BHS08213 | <b>BHS30 - 2100 (WP)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>PHOENIX multi-aesthetic</b>                                                      |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RP Ø 4,8              | BHS08411 | <b>BHS30 - 4100 (RP)</b>               | <br>     | <b>BHS03008</b>                    | <b>BHS03009</b>                    | <b>BHS02003</b> | BHS02004  | EM045169        |
| <b>GT MEDICAL COMPATIBLE</b>                                                        |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>BEST FIT HXE</b>                                                                 |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5              | BHS08211 | <b>BHS30 - 2100 (NP)</b>               | <br>     |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1              | BHS08212 | <b>BHS30 - 2100 (RP)</b>               | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1              | BHS08213 | <b>BHS30 - 2100 (WP)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>BEST FIT HXI</b>                                                                 |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | WP 4,5 / 5,0          | BHS08123 | <b>BHS30 - 1200 (WP Ø 4.5 / 5.0)</b>   |                                                                                                                                                                            |                                    |                                    | <b>BHS02009</b> | BHS02010  | EM046370        |
| <b>BEST FIT 002</b>                                                                 |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP Ø 3,5              | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
| <b>BEST FIT OCT</b>                                                                 |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RN Ø 4,8              | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b>         | <br> |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
|                                                                                     | SYN_RN Ø 4,8          | BHS08713 | <b>BHS30 - 7100 (RN SYNOCTA Ø 4,8)</b> | <br> |                                    | <b>BHS03019</b>                    | <b>BHS02013</b> | BHS02014  | EM046375        |
| <b>HI TECH COMPATIBLE</b>                                                           |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>SELF THREAD</b>                                                                  |                       |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | SLT Ø 3,3 / 3,75      | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | GTI Ø 4,2 / 5         | BHS08322 | <b>BHS30 - 3200 (RP)</b>               | <br> | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |

## 2. PRODUCTS

|                                                                                                                                                                         | PLATFORM Ø          | SCANBODY | LIBRARY PLATFORM               | <br> | CrCo BASE | Ti BASE  | CLINICAL SCREW | LAB SCREW | IMPLANT REPLICA |  |  |  |  |  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|----------|----------------|-----------|-----------------|--|--|--|--|--|
| LGI PLUS                                                                                                                                                                |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
|       | Ø 3,5 / 4,3 / 5 / 6 | BHS08321 | BHS30 - 3200 (NP)              |                                                                                       | BHS03050  | BHS03048 | BHS02021       | BHS02022  | EM045130        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                       | BHS03051  | BHS03049 |                |           |                 |  |  |  |  |  |
| IMPLANT DIRECT COMPATIBLE                                                                                                                                               |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
| INTERACTIVE                                                                                                                                                             |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
|       | NP Ø 3,5            | BHS08311 | BHS30 - 3100 (NP)              |                                                                                       | BHS03058  | BHS03056 | BHS02025       | BHS02026  | EM048148        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                       | BHS03059  | BHS03057 |                |           |                 |  |  |  |  |  |
|                                                                                                                                                                         | RP Ø 4,3            | BHS08312 | BHS30 - 3100 (RP)              |                                                                                       | BHS03046  | BHS03044 | BHS02023       | BHS02024  | EM048149        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                       | BHS03047  | BHS03045 |                |           |                 |  |  |  |  |  |
| LEGACY                                                                                                                                                                  |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
|       | NP Ø 3,5            | BHS08321 | BHS30 - 3200 (NP)              |                                                                                       | BHS03050  | BHS03048 | BHS02021       | BHS02022  | EM045130        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                       | BHS03051  | BHS03049 |                |           |                 |  |  |  |  |  |
|                                                                                                                                                                         | RP Ø 4,5            | BHS08322 | BHS30 - 3200 (RP)              |                                                                                       | BHS03033  | BHS03020 | BHS02021       | BHS02022  | EM046218        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                       | BHS03034  | BHS03021 |                |           |                 |  |  |  |  |  |
|                                                                                                                                                                         | WP Ø 5,7            | BHS08323 | BHS30 - 3200 (WP)              |                                                                                                                                                                        |           |          | BHS02021       | BHS02022  |                 |  |  |  |  |  |
| REPLUS / REPLANT                                                                                                                                                        |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
|   | NP Ø 3,5            | BHS08111 | BHS30 - 1100 (NP)              |                                                                                     |           | BHS03060 | BHS02027       | BHS02028  | EM046291        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                     |           | BHS03061 |                |           |                 |  |  |  |  |  |
|                                                                                                                                                                         | RP Ø 4,3            | BHS08112 | BHS30 - 1100 (RP)              |                                                                                     |           | BHS03012 | BHS02009       | BHS02010  | EM046292        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                     |           | BHS03013 |                |           |                 |  |  |  |  |  |
|                                                                                                                                                                         | WP Ø 5,0            | BHS08113 | BHS30 - 1100 (WP)              |                                                                                                                                                                        |           |          | BHS02009       | BHS02010  | EM046293        |  |  |  |  |  |
| SWISHPLUS / SWISHPLANT                                                                                                                                                  |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
|   | RN Ø 4,8            | BHS08712 | BHS30 - 7100 (RN Ø 4,8)        |                                                                                     |           | BHS03016 | BHS02011       | BHS02012  | EM045144        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                     |           | BHS03017 |                |           |                 |  |  |  |  |  |
|   | SYN_RN Ø 4,8        | BHS08713 | BHS30 - 7100 (RN SYNCTA Ø 4,8) |                                                                                     |           |          | BHS02013       | BHS02014  | EM046375        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                                                                                                        |           | BHS03019 |                |           |                 |  |  |  |  |  |
| JDEENTALCARE COMPATIBLE                                                                                                                                                 |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
| JDEVOLUTION                                                                                                                                                             |                     |          |                                |                                                                                                                                                                        |           |          |                |           |                 |  |  |  |  |  |
|   | Ø 3,7 / 4,3 / 5 / 6 | BHS08321 | BHS30 - 3200 (NP)              |                                                                                     | BHS03050  | BHS03048 | BHS02021       | BHS02022  | EM045130        |  |  |  |  |  |
|                                                                                                                                                                         |                     |          |                                |                                                                                     | BHS03051  | BHS03049 |                |           |                 |  |  |  |  |  |


## 2. PRODUCTS

|                                                                                     | PLATFORM Ø      | SCANBODY | LIBRARY PLATFORM                       |                                                                                           | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|-----------------|----------|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>KEYSTONE DENTAL COMPATIBLE</b>                                                   |                 |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>RESTORE</b>                                                                      |                 |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,3        | BHS08741 | <b>BHS30 - 7400 (S2)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02035</b> | BHS02036  | EM048022        |
|                                                                                     | RP Ø 4,1        | BHS08742 | <b>BHS30 - 7400 (S3 y S4)</b>          | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
| <b>OCT-XT</b>                                                                       |                 |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RN Ø 4,8        | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b>         | <br>     |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
|                                                                                     | SYN_RN Ø 4,8    | BHS08713 | <b>BHS30 - 7100 (RN SYNOCTA Ø 4,8)</b> |                                                                                           |                                    | <b>BHS03019</b>                    | <b>BHS02013</b> | BHS02014  | EM046375        |
| <b>RENOVA</b>                                                                       |                 |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,75       | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | RP Ø 4,5 / 4,75 | BHS08322 | <b>BHS30 - 3200 (RP)</b>               | <br>     | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
| <b>KLOCKNER COMPATIBLE</b>                                                          |                 |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>SK2 (EH)</b>                                                                     |                 |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | RP Ø 4,2        | BHS08521 | <b>BHS30 - 5200 (RP)</b>               | <br> |                                    |                                    | <b>BHS02031</b> | BHS02032  | EM046747        |
| <b>KL (EH)</b>                                                                      |                 |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP Ø 3,5        | BHS08211 | <b>BHS30 - 2100 (NP)</b>               | <br> |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1        | BHS08212 | <b>BHS30 - 2100 (RP)</b>               | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1        | BHS08213 | <b>BHS30 - 2100 (WP)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø               | SCANBODY | LIBRARY PLATFORM                             | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|--------------------------|----------|----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>LEADER COMPATIBLE</b>                                                            |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>TIXOS INTERNAL HEX</b>                                                           |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5 / 4,0           | BHS08321 | <b>BHS30 - 3200 (NP)</b>                     | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>TIXOS EXTERNAL HEX</b>                                                           |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RP Ø 4,1                 | BHS08612 | <b>BHS30 - 6100 (RP)</b>                     | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0                 | BHS08613 | <b>BHS30 - 6100 (WP)</b>                     |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>MEDENTIS COMPATIBLE</b>                                                          |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>ICX</b>                                                                          |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | Ø 3,45 / 3,7 - 4,1 / 4,8 | BHS08751 | <b>BHS30-7700 (Ø 3,45 / 3,7 - 4,1 / 4,8)</b> | <br>    | <b>BHS03180</b><br><b>BHS03181</b> | <b>BHS03178</b><br><b>BHS03179</b> | <b>BHS02015</b> | BHS02016  |                 |
|                                                                                     |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>MEGAGEN COMPATIBLE</b>                                                           |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>ANYONE INTERNAL</b>                                                              |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP                       | BHS08432 | <b>BHS30 - 4600 (4.0)</b>                    | <br> | <b>BHS03102</b><br><b>BHS03103</b> | <b>BHS03100</b>                    | <b>BHS02009</b> | BHS02010  |                 |
|                                                                                     | RP                       | BHS08432 | <b>BHS30 - 4600 (4.0)</b>                    | <br> | <b>BHS03102</b><br><b>BHS03103</b> | <b>BHS03100</b>                    | <b>BHS02009</b> | BHS02010  |                 |
|                                                                                     | WP                       | BHS08432 | <b>BHS30 - 4600 (4.0)</b>                    | <br> | <b>BHS03102</b><br><b>BHS03103</b> | <b>BHS03100</b>                    | <b>BHS02009</b> | BHS02010  |                 |
| <b>ANYONE EXTERNAL</b>                                                              |                          |          |                                              |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP                       | BHS08211 | <b>BHS30 - 2100 (NP)</b>                     | <br> |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP                       | BHS08212 | <b>BHS30 - 2100 (RP)</b>                     | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP                       | BHS08613 | <b>BHS30 - 6100 (WP)</b>                     |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø       | SCANBODY | LIBRARY PLATFORM                    | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|------------------|----------|-------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>ANYRIDGE</b>                                                                     |                  |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | Small Ø3,5       | BHS08991 | <b>BHS30-9900 (4,5/5,0)</b>         | <br>     | <b>BHS03166</b><br><b>BHS03167</b> | <b>BHS03164</b><br><b>BHS03165</b> | <b>BHS02027</b> | BHS02028  |                 |
|                                                                                     | Regular Ø4,0/4,5 | BHS08991 | <b>BHS30-9900 (4,5/5,0)</b>         | <br>     | <b>BHS03166</b><br><b>BHS03167</b> | <b>BHS03164</b><br><b>BHS03165</b> | <b>BHS02027</b> | BHS02028  |                 |
|                                                                                     | Wide Ø5,0/5,5    | BHS08991 | <b>BHS30-9900 (4,5/5,0)</b>         | <br>     | <b>BHS03166</b><br><b>BHS03167</b> | <b>BHS03164</b><br><b>BHS03165</b> | <b>BHS02027</b> | BHS02028  |                 |
| <b>OCTA Abutment</b>                                                                |                  |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | SYN_RN Ø4,8      | BHS08713 | <b>BHS30-7100 (RN SYNOCTA Ø4,8)</b> |                                                                                           |                                    | <b>BHS03019</b>                    | <b>BHS02013</b> | BHS02014  | EM046375        |
| <b>MICRODENT COMPATIBLE</b>                                                         |                  |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EKTOS (IH)</b>                                                                   |                  |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | NP Ø 3,5         | BHS08321 | <b>BHS30 - 3200 (NP)</b>            | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | RP Ø 4,5         | BHS08322 | <b>BHS30 - 3200 (RP)</b>            | <br>    | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
|                                                                                     | WP Ø 5,7         | BHS08323 | <b>BHS30 - 3200 (WP)</b>            |                                                                                                                                                                            |                                    |                                    | <b>BHS02021</b> | BHS02022  |                 |
| <b>MINI (EH)</b>                                                                    |                  |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP Ø 3,5         | BHS08611 | <b>BHS30 - 6100 (NP)</b>            |                                                                                                                                                                            |                                    |                                    | <b>BHS02031</b> | BHS02032  |                 |
| <b>SUV (EH)</b>                                                                     |                  |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | SUV Ø 4,1        | BHS08612 | <b>BHS30 - 6100 (RP)</b>            | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | SUV Ø 5,1        | BHS08613 | <b>BHS30 - 6100 (WP)</b>            |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>SYSTEM (EH)</b>                                                                  |                  |          |                                     |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | MST Ø 3,5        | BHS08611 | <b>BHS30 - 6100 (NP)</b>            |                                                                                                                                                                            |                                    |                                    | <b>BHS02031</b> | BHS02032  |                 |
|                                                                                     | MS Ø 3,5         | BHS08611 | <b>BHS30 - 6100 (NP)</b>            |                                                                                                                                                                            |                                    |                                    | <b>BHS02031</b> | BHS02032  |                 |



## 2. PRODUCTS

|                                                                                     | PLATFORM Ø            | SCANBODY | LIBRARY PLATFORM                     | <br>     | CrCo<br>BASE                       | Ti<br>BASE                         | CLINICAL<br>SCREW | LAB<br>SCREW | IMPLANT<br>REPLICA |
|-------------------------------------------------------------------------------------|-----------------------|----------|--------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-------------------|--------------|--------------------|
| <b>UNIVERSAL (EH)</b>                                                               |                       |          |                                      |                                                                                                                                                                            |                                    |                                    |                   |              |                    |
|    | MST Ø 3,5             | BHS08611 | <b>BHS30 - 6100 (NP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02031</b>   | BHS02032     |                    |
|                                                                                     | MS Ø 3,5              | BHS08611 | <b>BHS30 - 6100 (NP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02031</b>   | BHS02032     |                    |
|                                                                                     | SU Ø 4,1              | BHS08612 | <b>BHS30 - 6100 (RP)</b>             | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b>   | BHS02002     | EM045098           |
|                                                                                     | SU Ø 5,1              | BHS08613 | <b>BHS30 - 6100 (WP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b>   | BHS02002     | EM045082           |
| <b>MIS COMPATIBLE</b>                                                               |                       |          |                                      |                                                                                                                                                                            |                                    |                                    |                   |              |                    |
| <b>LANCE</b>                                                                        |                       |          |                                      |                                                                                                                                                                            |                                    |                                    |                   |              |                    |
|    | STD Ø 4,1             | BHS08612 | <b>BHS30 - 6100 (RP)</b>             | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b>   | BHS02002     | EM045098           |
|                                                                                     | WP Ø 5,0              | BHS08613 | <b>BHS30 - 6100 (WP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b>   | BHS02002     | EM045082           |
| <b>SEVEN / M4</b>                                                                   |                       |          |                                      |                                                                                                                                                                            |                                    |                                    |                   |              |                    |
|   | Narrow Ø 3,30         | BHS08951 | <b>BHS30-9500 (NARROW)</b>           | <br>    | <b>BHS03200</b><br><b>BHS03201</b> | <b>BHS03198</b><br><b>BHS03199</b> | <b>BHS02041</b>   | BHS02042     |                    |
|                                                                                     | Standard Ø 3,75 / 4,2 | BHS08321 | <b>BHS30 - 3200 (NP)</b>             | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b>   | BHS02022     | EM045130           |
|                                                                                     | Wide Ø 5,0 / 6,0      | BHS08322 | <b>BHS30 - 3200 (RP)</b>             | <br> | <b>BHS03033</b>                    | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b>   | BHS02022     | EM046218           |
| <b>MULTIUNIT MIS</b>                                                                |                       |          |                                      |                                                                                                                                                                            |                                    |                                    |                   |              |                    |
|  | RP Ø 4,8              | BHS08411 | <b>BHS30 - 4100 (RP)</b>             |                                                                                         | <b>BHS03008</b>                    | <b>BHS03009</b>                    | <b>BHS02037</b>   | BHS02038     |                    |
| <b>MOZO GRAU COMPATIBLE</b>                                                         |                       |          |                                      |                                                                                                                                                                            |                                    |                                    |                   |              |                    |
| <b>IN HEX</b>                                                                       |                       |          |                                      |                                                                                                                                                                            |                                    |                                    |                   |              |                    |
|  | RP Ø 3,75 / 4,25      | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 / 4.0)</b> | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b>   | BHS02026     | EM046371           |
|                                                                                     | WP Ø 5,0              | BHS08123 | <b>BHS30 - 1200 (WP Ø 4.5 / 5.0)</b> |                                                                                                                                                                            |                                    |                                    | <b>BHS02009</b>   | BHS02010     | EM046370           |
|                                                                                     | CONO 20               | BHS08124 | <b>BHS30 - 1200 (20°)</b>            |                                                                                         |                                    | <b>BHS03072</b>                    | <b>BHS02003</b>   | BHS02004     | EM045638           |
|                                                                                     | CONO 45               | BHS08125 | <b>BHS30 - 1200 (45°)</b>            |                                                                                         |                                    | <b>BHS03074</b>                    | <b>BHS02003</b>   | BHS02004     | EM047786           |


## 2. PRODUCTS

|                                                                                     | PLATFORM Ø          | SCANBODY | LIBRARY PLATFORM                     | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|---------------------|----------|--------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>OSSEUS</b>                                                                       |                     |          |                                      |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,4            | BHS08611 | <b>BHS30 - 6100 (NP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045468        |
|                                                                                     | RP Ø 4,1            | BHS08612 | <b>BHS30 - 6100 (RP)</b>             | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0            | BHS08213 | <b>BHS30 - 2100 (WP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>MPI COMPATIBLE</b>                                                               |                     |          |                                      |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EXTERNA EH</b>                                                                   |                     |          |                                      |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5            | BHS08421 | <b>BHS30 - 4200 (NP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM047957        |
|                                                                                     | RP Ø 4,1            | BHS08612 | <b>BHS30 - 6100 (RP)</b>             | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0            | BHS08613 | <b>BHS30 - 6100 (WP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>INTERNA CM</b>                                                                   |                     |          |                                      |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|   | RP Ø 4,0            | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 / 4.0)</b> | <br>    | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
|                                                                                     | WP Ø 5,0            | BHS08123 | <b>BHS30 - 1200 (WP Ø 4.5 / 5.0)</b> |                                                                                                                                                                            |                                    |                                    | <b>BHS02009</b> | BHS02010  | EM046370        |
| <b>NEOBIOTECH COMPATIBLE</b>                                                        |                     |          |                                      |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>EB EXTERNAL SYSTEM</b>                                                           |                     |          |                                      |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | NP Ø 3,5            | BHS08211 | <b>BHS30 - 2100 (NP)</b>             | <br> |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1            | BHS08212 | <b>BHS30 - 2100 (RP)</b>             | <br> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1            | BHS08213 | <b>BHS30 - 2100 (WP)</b>             |                                                                                                                                                                            |                                    |                                    | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>IS INTERNAL SYSTEM</b>                                                           |                     |          |                                      |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 3,5 / 4 / 4,5 / 5 | BHS08432 | <b>BHS30 - 4600 (4.0)</b>            | <br> | <b>BHS03102</b><br><b>BHS03103</b> | <b>BHS03100</b><br><b>BHS03006</b> | <b>BHS02009</b> | BHS02010  |                 |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø       | SCANBODY | LIBRARY PLATFORM         | <div> <div>○</div> <div>○</div> </div> | CrCo BASE                                            | Ti BASE                                              | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|------------------|----------|--------------------------|----------------------------------------|------------------------------------------------------|------------------------------------------------------|-----------------|-----------|-----------------|
| <b>NEODENT COMPATIBLE</b>                                                           |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
| <b>TITAMAX HE SMART</b>                                                             |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
|    | RP Ø 4,1         | BHS08612 | <b>BHS30 - 6100 (RP)</b> | <div> <div>○</div> <div>○</div> </div> | <div> <div>BHS03003</div> <div>BHS03004</div> </div> | <div> <div>BHS03005</div> <div>BHS03006</div> </div> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0         | BHS08613 | <b>BHS30 - 6100 (WP)</b> |                                        |                                                      |                                                      | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>MINI PILAR CM</b>                                                                |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
|    | RP Ø 4,8         | BHS08411 | <b>BHS30 - 4100 (RP)</b> | <div> <div>○</div> </div>              | <b>BHS03008</b>                                      | <b>BHS03009</b>                                      | <b>BHS02003</b> | BHS02004  | EM045169        |
| <b>NOBEL COMPATIBLE</b>                                                             |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
| <b>BRÄNEMARK</b>                                                                    |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
|   | NP Ø 3,5         | BHS08211 | <b>BHS30 - 2100 (NP)</b> | <div> <div>○</div> <div>○</div> </div> |                                                      | <div> <div>BHS03010</div> <div>BHS03011</div> </div> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1         | BHS08212 | <b>BHS30 - 2100 (RP)</b> | <div> <div>○</div> <div>○</div> </div> | <div> <div>BHS03003</div> <div>BHS03004</div> </div> | <div> <div>BHS03005</div> <div>BHS03006</div> </div> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,1         | BHS08213 | <b>BHS30 - 2100 (WP)</b> |                                        |                                                      |                                                      | <b>BHS02039</b> | BHS02040  | EM045070        |
| <b>REPLACE</b>                                                                      |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
|  | NP Ø 3,5         | BHS08111 | <b>BHS30 - 1100 (NP)</b> | <div> <div>○</div> <div>○</div> </div> |                                                      | <div> <div>BHS03060</div> <div>BHS03061</div> </div> | <b>BHS02027</b> | BHS02028  | EM046291        |
|                                                                                     | RP Ø 4,3         | BHS08112 | <b>BHS30 - 1100 (RP)</b> | <div> <div>○</div> <div>○</div> </div> |                                                      | <div> <div>BHS03012</div> <div>BHS03013</div> </div> | <b>BHS02009</b> | BHS02010  | EM046292        |
|                                                                                     | WP Ø 5,0         | BHS08113 | <b>BHS30 - 1100 (WP)</b> |                                        |                                                      |                                                      | <b>BHS02009</b> | BHS02010  | EM046293        |
| <b>ACTIVE</b>                                                                       |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
|  | NP Ø 3,5         | BHS08311 | <b>BHS30 - 3100 (NP)</b> | <div> <div>○</div> <div>○</div> </div> | <div> <div>BHS03058</div> <div>BHS03059</div> </div> | <div> <div>BHS03056</div> <div>BHS03057</div> </div> | <b>BHS02025</b> | BHS02026  | EM048148        |
|                                                                                     | RP Ø 4,3 / Ø 5,0 | BHS08312 | <b>BHS30 - 3100 (RP)</b> | <div> <div>○</div> <div>○</div> </div> | <div> <div>BHS03046</div> <div>BHS03047</div> </div> | <div> <div>BHS03044</div> <div>BHS03045</div> </div> | <b>BHS02023</b> | BHS02024  | EM048149        |
| <b>MULTIUNIT</b>                                                                    |                  |          |                          |                                        |                                                      |                                                      |                 |           |                 |
|  | RP Ø 4,8         | BHS08411 | <b>BHS30 - 4100 (RP)</b> | <div> <div>○</div> </div>              | <b>BHS03008</b>                                      | <b>BHS03009</b>                                      | <b>BHS02003</b> | BHS02004  | EM045169        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø        |       | SCANBODY | LIBRARY PLATFORM   |                                                                                           | CrCo BASE            | Ti BASE              | CLINICAL SCREW | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|-------------------|-------|----------|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|----------------------|----------------|-----------|-----------------|
| NT TRADING COMPATIBLE                                                               |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
| 2-CONnect                                                                           |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
|    | M1,4              | Ø5,0  | BHS08771 | BHS30-7700 (M1,4)  |                                                                                           | BHS03177             | BHS03176             | BHS02063       | BHS02064  |                 |
|                                                                                     | M1,6              | Ø5,0  | BHS08771 | BHS30-7700 (M1,6)  |                                                                                           | BHS03177             | BHS03176             | BHS02065       | BHS02066  |                 |
|                                                                                     | M1,8              | Ø5,0  | BHS08771 | BHS30-7700 (M1,8)  |                                                                                           | BHS03177             | BHS03176             | BHS02067       | BHS02068  |                 |
| OSSTEM IMPLANT COMPATIBLE                                                           |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
| OSSTEM TS                                                                           |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
|    | Ø 4 / 4,5 / 5 / 6 |       | BHS08432 | BHS30 - 4600 (4.0) | <br>     | BHS03102<br>BHS03103 | BHS03100<br>BHS03101 | BHS02009       | BHS02010  |                 |
| OSSTEM US                                                                           |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | MINI              | Ø 3,5 | BHS08211 | BHS30 - 2100 (NP)  | <br> |                      | BHS03010<br>BHS03011 | BHS02005       | BHS02006  | EM045222        |
|                                                                                     | RP                | Ø 4,1 | BHS08212 | BHS30 - 2100 (RP)  | <br> | BHS03003<br>BHS03004 | BHS03005<br>BHS03006 | BHS02001       | BHS02002  | EM045098        |
|                                                                                     | WP                | Ø 5,1 | BHS08213 | BHS30 - 2100 (WP)  |                                                                                                                                                                            |                      |                      | BHS02039       | BHS02040  | EM045070        |
| OSTEOPLUS COMPATIBLE                                                                |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
| ARGENTIUM                                                                           |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | RP                | Ø 4,2 | BHS08521 | BHS30 - 5200 (RP)  |                                                                                                                                                                            |                      |                      | BHS02031       | BHS02032  | EM046747        |
| OSTEOPLUS - SHc                                                                     |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | RP                | Ø 4,1 | BHS08612 | BHS30 - 6100 (RP)  | <br> | BHS03003<br>BHS03004 | BHS03005<br>BHS03006 | BHS02001       | BHS02002  | EM045098        |
|                                                                                     | WP                | Ø 5,0 | BHS08613 | BHS30 - 6100 (WP)  |                                                                                                                                                                            |                      |                      | BHS02001       | BHS02002  | EM045082        |
| OSTEOPLUS - SHi                                                                     |                   |       |          |                    |                                                                                                                                                                            |                      |                      |                |           |                 |
|  | RP                | Ø 3,5 | BHS08321 | BHS30 - 3200 (NP)  | <br> | BHS03050<br>BHS03051 | BHS03048<br>BHS03049 | BHS02021       | BHS02022  | EM045130        |
|                                                                                     | WP                | Ø 4,5 | BHS08322 | BHS30 - 3200 (RP)  | <br> | BHS03033<br>BHS03034 | BHS03020<br>BHS03021 | BHS02021       | BHS02022  | EM046218        |



## 2. PRODUCTS

|                                                                                     | PLATFORM Ø                  | SCANBODY | LIBRARY PLATFORM                       | <br>     | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|-----------------------------|----------|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>OXTEIN COMPATIBLE</b>                                                            |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>HEX EXTERNA</b>                                                                  |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | MINI Ø 3,3                  | BHS08741 | <b>BHS30 - 7400 (S2)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02035</b> | BHS02036  | EM048022        |
|                                                                                     | RP Ø 4,1                    | BHS08742 | <b>BHS30 - 7400 (S3 y S4)</b>          | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | WP Ø 5,0                    | BHS08613 | <b>BHS30 - 6100 (WP)</b>               |                                                                                                                                                                            |                                    |                                    | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>HEX INTERNA</b>                                                                  |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5                    | BHS08321 | <b>BHS30 - 3200 (NP)</b>               | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | RP Ø 4,5                    | BHS08322 | <b>BHS30 - 3200 (RP)</b>               | <br>     | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
| <b>OCTA INTERNA</b>                                                                 |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|    | RN Ø 4,8                    | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b>         | <br>     |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
|                                                                                     | RN Ø 4,8                    | BHS08713 | <b>BHS30 - 7100 (RN SYNOCTA Ø 4,8)</b> |                                                                                           |                                    | <b>BHS03019</b>                    | <b>BHS02013</b> | BHS02014  | EM046375        |
| <b>PALTOP COMPATIBLE</b>                                                            |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
| <b>INTERNAL HEX CONNECTION</b>                                                      |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Narrow Ø 3,25               | BHS08951 | <b>BHS30-9500 (NARROW)</b>             | <br> | <b>BHS03200</b><br><b>BHS03201</b> | <b>BHS03198</b><br><b>BHS03199</b> | <b>BHS02041</b> | BHS02042  |                 |
|                                                                                     | Standard Ø 3,75 / 4,2 / 5,0 | BHS08321 | <b>BHS30-3200 (NP)</b>                 | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | Wide Ø 6,0                  | BHS08322 | <b>BHS30-3200 (RP)</b>                 | <br> | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
| <b>PCIA Implant System (Paltop Conical Active)</b>                                  |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Narrow Ø 3,25               | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 /4.0)</b>    | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
|                                                                                     | Standard Ø 3,75 / 4,2 / 5,0 | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 /4.0)</b>    | <br> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
| <b>MULTI UNIT</b>                                                                   |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 5,0                       | BHS08271 | <b>BHS30-2700 (MU)</b>                 |                                                                                         |                                    | <b>BHS03194</b>                    | <b>BHS02031</b> | BHS02032  |                 |
| <b>SINGLE UNIT</b>                                                                  |                             |          |                                        |                                                                                                                                                                            |                                    |                                    |                 |           |                 |
|  | Ø 4,5                       | BHS08251 | <b>BHS30-2500 (SINGLE)</b>             |                                                                                         |                                    | <b>BHS03196</b>                    | <b>BHS02031</b> | BHS02032  |                 |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø    | SCANBODY | LIBRARY PLATFORM                     | <div> <div>○</div> <div>○</div> </div> | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|---------------|----------|--------------------------------------|----------------------------------------|------------------------------------|------------------------------------|-----------------|-----------|-----------------|
| <b>PHIBO COMPATIBLE</b>                                                             |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
| TSH                                                                                 |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
|    | S2 Ø 3,3      | BHS08741 | <b>BHS30 - 7400 (S2)</b>             |                                        |                                    |                                    | <b>BHS02035</b> | BHS02036  | EM048022        |
|                                                                                     | S3 / S4 Ø 4,1 | BHS08742 | <b>BHS30 - 7400 (S3 y S4)</b>        | <div> <div>○</div> <div>○</div> </div> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
| AUREA EVO                                                                           |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5      |          |                                      | <div> <div>○</div> <div>○</div> </div> | <b>BHS03106</b><br><b>BHS03107</b> |                                    | <b>BHS02041</b> | BHS02042  |                 |
|                                                                                     | RP Ø 4,3      |          |                                      | <div> <div>○</div> <div>○</div> </div> | <b>BHS03110</b><br><b>BHS03111</b> |                                    | <b>BHS02021</b> | BHS02022  |                 |
| <b>RADHEX IMPLANTS COMPATIBLE</b>                                                   |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
| PHE / PHEA                                                                          |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
|    | NP Ø 3,5      | BHS08211 | <b>BHS30 - 2100 (NP)</b>             | <div> <div>○</div> <div>○</div> </div> |                                    | <b>BHS03010</b><br><b>BHS03011</b> | <b>BHS02005</b> | BHS02006  | EM045222        |
|                                                                                     | RP Ø 4,1      | BHS08212 | <b>BHS30 - 2100 (RP)</b>             | <div> <div>○</div> <div>○</div> </div> | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
| PHI / PHIA                                                                          |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
|  | NP Ø3,5       | BHS08321 | <b>BHS30-3200 (NP)</b>               | <div> <div>○</div> <div>○</div> </div> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     | RP Ø4,5       | BHS08322 | <b>BHS30-3200 (RP)</b>               | <div> <div>○</div> <div>○</div> </div> | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | BHS02022  | EM046218        |
| PCI                                                                                 |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
|  | 280           | BHS08122 | <b>BHS30 - 1200 (RP Ø 3.5 / 4.0)</b> | <div> <div>○</div> <div>○</div> </div> | <b>BHS03070</b><br><b>BHS03071</b> | <b>BHS03068</b><br><b>BHS03069</b> | <b>BHS02025</b> | BHS02026  | EM046371        |
|                                                                                     | 350           | BHS08123 | <b>BHS30 - 1200 (WP Ø 4.5 / 5.0)</b> |                                        |                                    |                                    | <b>BHS02009</b> | BHS02010  | EM046370        |
| <b>SGS COMPATIBLE</b>                                                               |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
| P1 / P7                                                                             |               |          |                                      |                                        |                                    |                                    |                 |           |                 |
|  | RP Ø 3,75     | BHS08321 | <b>BHS30 - 3200 (NP)</b>             | <div> <div>○</div> <div>○</div> </div> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | BHS02022  | EM045130        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø         | SCANBODY | LIBRARY PLATFORM                          | <br>CrCo<br>BASE                                                                          | Ti<br>BASE                         | CLINICAL<br>SCREW                  | LAB<br>SCREW    | IMPLANT<br>REPLICA   |
|-------------------------------------------------------------------------------------|--------------------|----------|-------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|----------------------|
| <b>SIGNO VINCES COMPATIBLE</b>                                                      |                    |          |                                           |                                                                                                                                                                            |                                    |                                    |                 |                      |
| DUO / INTEGRA                                                                       |                    |          |                                           |                                                                                                                                                                            |                                    |                                    |                 |                      |
|    | RP Ø 4,1           | BHS08212 | <b>BHS30-2100</b><br>(RP)                 | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | BHS02002<br>EM045098 |
| DUOCON / CM                                                                         |                    |          |                                           |                                                                                                                                                                            |                                    |                                    |                 |                      |
|    | WP_4.5 / 5.0       | BHS08123 | <b>BHS30-1200</b><br>(WP Ø 4.5 / 5.0)     |                                                                                                                                                                            |                                    | <b>BHS02009</b>                    | BHS02010        | EM046370             |
| <b>STRAUMANN COMPATIBLE</b>                                                         |                    |          |                                           |                                                                                                                                                                            |                                    |                                    |                 |                      |
| OCTA INTERN                                                                         |                    |          |                                           |                                                                                                                                                                            |                                    |                                    |                 |                      |
|    | RN Ø 4,8           | BHS08712 | <b>BHS30 - 7100</b><br>(RN Ø 4,8)         | <br>     | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b>                    | BHS02012        | EM045144             |
|  | SYN_RN Ø 4,8       | BHS08713 | <b>BHS30 - 7100</b><br>(RN SYNOCTA Ø 4,8) |                                                                                          | <b>BHS03019</b>                    | <b>BHS02013</b>                    | BHS02014        | EM046375             |
| BONE LEVEL                                                                          |                    |          |                                           |                                                                                                                                                                            |                                    |                                    |                 |                      |
|  | NC Ø 3,3           | BHS08811 | <b>BHS30 - 8100</b><br>(NC)               | <br> | <b>BHS03064</b><br><b>BHS03065</b> | <b>BHS02015</b>                    | BHS02016        | EM048048             |
|                                                                                     | RC Ø 4,1/<br>Ø 4,8 | BHS08812 | <b>BHS30 - 8100</b><br>(RC)               | <br> | <b>BHS03036</b><br><b>BHS03037</b> | <b>BHS02015</b>                    | BHS02016        | EM046904             |
| BL MULTIBASE                                                                        |                    |          |                                           |                                                                                                                                                                            |                                    |                                    |                 |                      |
|  | RC Ø 4,8           |          |                                           |                                                                                                                                                                            |                                    | <b>BHS02017</b>                    | BHS02018        | EM047780             |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø         | SCANBODY | LIBRARY PLATFORM                        |    | CrCo BASE       | Ti BASE         | CLINICAL SCREW  | LAB SCREW | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|--------------------|----------|-----------------------------------------|-------------------------------------------------------------------------------------|-----------------|-----------------|-----------------|-----------|-----------------|
| <b>SWEDEN &amp; MARTINA COMPATIBLE</b>                                              |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
| <b>PREMIUM KOHNO</b>                                                                |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
|    | NP Ø 3,3           | BHS08961 | <b>BHS30 - 9600 (NP Ø 3,3)</b>          |                                                                                     |                 |                 | <b>BHS02021</b> | BHS02022  |                 |
|                                                                                     | RP Ø 3,8           | BHS08962 | <b>BHS30 - 9600 (RP Ø 3,8)</b>          |                                                                                     |                 |                 | <b>BHS02021</b> | BHS02022  |                 |
|                                                                                     | WP Ø 4,25          | BHS08963 | <b>BHS30 - 9600 (WP Ø 4,25)</b>         |                                                                                     |                 |                 | <b>BHS02011</b> | BHS02012  |                 |
| <b>SYRA (EH)</b>                                                                    |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
|    | Ø 4,10 / 4,25      | BHS08612 | <b>BHS30 - 6100 (RP)</b>                |    | <b>BHS03003</b> | <b>BHS03005</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | Ø 5,0              | BHS08613 | <b>BHS30 - 6100 (WP)</b>                |    | <b>BHS03004</b> | <b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>SYRA Short (EH)</b>                                                              |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
|    | Ø 4,10             | BHS08612 | <b>BHS30 - 6100 (RP)</b>                |    | <b>BHS03003</b> | <b>BHS03005</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     | Ø 5,0              | BHS08613 | <b>BHS30 - 6100 (WP)</b>                |    | <b>BHS03004</b> | <b>BHS03006</b> | <b>BHS02001</b> | BHS02002  | EM045082        |
| <b>TRI DENTAL IMPLANTS COMPATIBLE</b>                                               |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
| <b>TRI - VENT</b>                                                                   |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
|  | Ø 3,75 / 4,1 / 4,7 | BHS08321 | <b>BHS30 - 3200 (NP)</b>                |  | <b>BHS03050</b> | <b>BHS03048</b> | <b>BHS02021</b> | BHS02022  | EM045130        |
|                                                                                     |                    |          |                                         |  | <b>BHS03051</b> | <b>BHS03049</b> |                 |           |                 |
| <b>TRINON COMPATIBLE</b>                                                            |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
| <b>TRINON - Q2</b>                                                                  |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
|  | RP Ø 4,1           | BHS08212 | <b>BHS30 - 2100 (RP)</b>                |  | <b>BHS03003</b> | <b>BHS03005</b> | <b>BHS02001</b> | BHS02002  | EM045098        |
|                                                                                     |                    |          |                                         |  | <b>BHS03004</b> | <b>BHS03006</b> |                 |           |                 |
| <b>TRINON - QK</b>                                                                  |                    |          |                                         |                                                                                     |                 |                 |                 |           |                 |
|  | RN Ø 4,8           | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b>          |  |                 | <b>BHS03016</b> | <b>BHS02011</b> | BHS02012  | EM045144        |
|                                                                                     |                    |          |                                         |  |                 | <b>BHS03017</b> |                 |           |                 |
|                                                                                     | RN Ø 4,8           | BHS08713 | <b>BHS30 - 7100 (RN SYNOC TA Ø 4,8)</b> |  |                 | <b>BHS03019</b> | <b>BHS02013</b> | BHS02014  | EM046375        |

## 2. PRODUCTS

|                                                                                     | PLATFORM Ø | SCANBODY | LIBRARY PLATFORM                        |                                                                                           | CrCo BASE                          | Ti BASE                            | CLINICAL SCREW  | LAB SCREW       | IMPLANT REPLICA |
|-------------------------------------------------------------------------------------|------------|----------|-----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|-----------------|-----------------|-----------------|
| <b>ZIACOM COMPATIBLE</b>                                                            |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
| <b>ZINIC - IH</b>                                                                   |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
|    | NP Ø 3,2   | BHS08951 | <b>BHS30-9500 (NARROW)</b>              | <br>     | <b>BHS03200</b><br><b>BHS03201</b> | <b>BHS03198</b><br><b>BHS03199</b> | <b>BHS02041</b> | <b>BHS02042</b> |                 |
|                                                                                     | RP Ø 3,5   | BHS08321 | <b>BHS30 - 3200 (NP)</b>                | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | <b>BHS02022</b> | EM045130        |
|                                                                                     | WP Ø 4,5   | BHS08322 | <b>BHS30 - 3200 (RP)</b>                | <br>     | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | <b>BHS02022</b> | EM046218        |
| <b>OEX - EH</b>                                                                     |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
|    | NP Ø 3,3   | BHS08741 | <b>BHS30 - 7400 (S2)</b>                |                                                                                                                                                                            |                                    |                                    | <b>BHS02035</b> | <b>BHS02036</b> | EM048022        |
|                                                                                     | RP Ø 4,1   | BHS08212 | <b>BHS30 - 2100 (RP)</b>                | <br>     | <b>BHS03003</b><br><b>BHS03004</b> | <b>BHS03005</b><br><b>BHS03006</b> | <b>BHS02001</b> | <b>BHS02002</b> | EM045098        |
| <b>OIN - IH (M2)</b>                                                                |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
|    | RP Ø 4,0   | BHS08321 | <b>BHS30-3200 (NP)</b>                  | <br>     | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02011</b> | <b>BHS02012</b> |                 |
|                                                                                     | WP Ø 5,0   | BHS08322 | <b>BHS30-3200 (RP)</b>                  | <br>     | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02011</b> | <b>BHS02012</b> |                 |
| <b>XDrive (Multi-unit)</b>                                                          |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
|   | RP Ø 4,8   | BHS08411 | <b>BHS30 - 4100 (RP)</b>                |                                                                                         | <b>BHS03008</b><br><b>BHS03009</b> |                                    | <b>BHS02003</b> | <b>BHS02004</b> | EM045169        |
| <b>ZIMMER COMPATIBLE</b>                                                            |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
| <b>INTERNAL HEX</b>                                                                 |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
|  | NP Ø 3,5   | BHS08321 | <b>BHS30 - 3200 (NP)</b>                | <br> | <b>BHS03050</b><br><b>BHS03051</b> | <b>BHS03048</b><br><b>BHS03049</b> | <b>BHS02021</b> | <b>BHS02022</b> | EM045130        |
|                                                                                     | RP Ø 4,5   | BHS08322 | <b>BHS30 - 3200 (RP)</b>                | <br> | <b>BHS03033</b><br><b>BHS03034</b> | <b>BHS03020</b><br><b>BHS03021</b> | <b>BHS02021</b> | <b>BHS02022</b> | EM046218        |
|                                                                                     | WP Ø 5,7   | BHS08323 | <b>BHS30 - 3200 (WP)</b>                |                                                                                                                                                                            |                                    |                                    | <b>BHS02021</b> | <b>BHS02022</b> |                 |
| <b>SWISS - PLUS</b>                                                                 |            |          |                                         |                                                                                                                                                                            |                                    |                                    |                 |                 |                 |
|  | RP Ø 4,8   | BHS08712 | <b>BHS30 - 7100 (RN Ø 4,8)</b>          | <br> |                                    | <b>BHS03016</b><br><b>BHS03017</b> | <b>BHS02011</b> | <b>BHS02012</b> | EM045144        |
|                                                                                     | RP Ø 4,8   | BHS08713 | <b>BHS30 - 7100 (RN SYNOC TA Ø 4,8)</b> |                                                                                         |                                    | <b>BHS03019</b>                    | <b>BHS02013</b> | <b>BHS02014</b> | EM046375        |

### 3. BHS30 SCREWDRIVER

---



#### EXCELLENT MECHANICAL PROPERTIES & SCREW HEAD PROTECTION

---

Universal contra angle compatible connection (ISO1797).

Patented design offering excellent torsion resistance (90 N-cm).

Optimal coupling engagement ensuring screw head integrity, guaranteed for life.

TiN Coated: Surface treatment which increases the life of the driver.

Designed to correct 30° angle deviations with 100% efficient torque transmission.

#### SCREWDRIVER BHS

---



BHS01001 - L25mm  
BHS01004 - L20mm  
BHS01005 - L32mm

#### HANDLER

---



BHS01110 - Nob Ø 8mm / CA  
BHS01111 - Str Ø 7mm / CA

## 4. PROTOCOLS


### CAD/CAM IMPLANT LEVEL

#### Products required:


BHS30 Screw  
BHS30 Screwdriver  
BHS30 Scanbody and associated BHS30 Cad-Cam Library  
BHS30 T-Tool

#### Materials:


Metal Sintering or Milling




**1.**  
CLINIC\_ Assembly of the BHS30 IO Scanbody on the patient's implants (flat surface aligned to vestibular) and Intraoral scanning. Delivery of the resulting .stl file to the laboratory.




**2.**  
LAB\_ Assembly of the BHS EO Scanbody on the implant analogues (flat surface aligned to vestibular) and model scanning in accordance with the software protocols to obtain the .stl file.




**3.**  
LAB\_ CAD: Design of the Implant Level structure from .stl file with the corresponding BHS30 Cad Library and selection of the desired Angulation -Free Angulation- for the BHS Screw access channels, according to the undesired emergence on the restoration.



**4.**  
LAB\_ CAM: Production of the structure as per the software protocols. For a perfect fit and precise finishing of the BHS30 Screw seating platform, the use of the BHS30 T-Tool is recommended. Please refer to our Milling Support service to check on software and T-Tool compatibilities.




**5.**  
LAB\_ Adjustment of the structure on the analogues until obtaining passive fit on the model.




**6.**  
LAB\_ Adjustment of the structure on the analogues until obtaining passive fit on the model.

## 4. PROTOCOLS


### CAD/CAM IMPLANT LEVEL



**7.**  
CLINIC\_ Metal test, especially in case of multiple structures, to confirm adjustment on patient's implants and their positioning in the mouth.



**8.**  
LAB\_ Final adjustments in the structure, where appropriate, aesthetic restoration, polishing and finishing.



**9.**  
CLINIC\_ Placement of the prosthesis to the patient.

## 4. PROTOCOLS


### CAD/CAM SINTERING TO CrCo BASE

#### Products required:


BHS30 Screw  
BHS30 Screwdriver  
BHS30 Scanbody and associated BHS30 Cad-Cam Library  
BHS30 CrCo Base

#### Materials:


CrCo Sintering or Milling




**1.**  
CLINIC\_ Assembly of the BHS30 IO Scanbody on the patient's implants (flat surface aligned to vestibular) and Intraoral scanning. Delivery of the resulting .stl file to the laboratory.




**2.**  
LAB\_ Assembly of the BHS EO Scanbody on the implant analogues (flat surface aligned to vestibular) and model scanning in accordance with the software protocols to obtain the .stl file.




**3.**  
LAB\_ CAD: Design of the structure onto BHS30 CrCo Base from .stl file with the corresponding BHS30 Cad Library and selection of the desired Angulation -Predefined Angulation- for the BHS30 Screw access channels, according to the undesired emergence on the restoration.



**4.**  
LAB\_ CAM: Production of the structure onto BHS30 CrCo Base as per the software protocols. Please refer to our Milling Support service to check on software compatibilities.




**5.**  
LAB\_ Adjustment of the metal cast structure on the BHS30 CrCo Base until obtaining passive fit on the model and perimetral laser welding of both parts.




**6.**  
LAB\_ Adjustment of the metal cast structure on the BHS30 CrCo Base until obtaining passive fit on the model and sending the Base and the structure to the Clinic separately for metal testing according to the actual relation in the mouth.

## 4. PROTOCOLS


### CAD/CAM SINTERING TO CrCo BASE



**7.**  
CLINIC\_ Metal test, especially in case of multiple structures, to confirm adjustment on patient's implants and their positioning in the mouth. Provisional fixation with passive fit in the mouth and delivery back to the lab.



**8.**  
LAB\_ Final adjustments in the structure, where appropriate, aesthetic restoration, polishing and finishing. Reassembly of the replicas according to the actual relation in the mouth and perimetral laser welding of both parts. Aesthetic restoration, polishing and finishing.



**9.**  
CLINIC\_ Placement of the prosthesis to the patient.

## 4. PROTOCOLS


### CAD/CAM MILLING TO BHS30 Ti BASE

#### Products required:


BHS30 Screw  
BHS30 Screwdriver  
BHS30 Scanbody and associated BHS30 Cad-Cam Library  
BHS30 CrCo Base

#### Materials:


Metal Sintering or Milling




**1.**  
CLINIC\_ Assembly of the BHS30 IO Scanbody on the patient's implants (flat surface aligned to vestibular) and Intraoral scanning. Delivery of the resulting .stl file to the laboratory.




**2.**  
LAB\_ Assembly of the BHS EO Scanbody on the implant analogues (flat surface aligned to vestibular) and model scanning in accordance with the software protocols to obtain the .stl file.




**3.**  
LAB\_ CAD: Design of the structure onto BHS30 Ti Base from .stl file with the corresponding BHS30 Cad Library and selection of the desired Angulation -Predefined Angulation- for the BHS30 Screw access channels, according to the undesired emergence on the restoration.



**4.**  
LAB\_ CAM: Production of the structure onto BHS30 CrCo Base as per the software protocols. Please refer to our Milling Support service to check on software compatibilities.




**5.**  
LAB\_ Adjustment of the metal cast structure on the BHS30 Ti Base until obtaining passive fit on the model and cementing of both parts.




**6.**  
LAB\_ Adjustment of the metal cast structure on the BHS30 CrCo Base until obtaining passive fit on the model and sending the Base and the structure to the Clinic separately for metal testing according to the actual relation in the mouth.

## 4. PROTOCOLS


### CAD/CAM MILLING TO Ti BASE



**7.**  
CLINIC\_ Restoration test, especially in case of multiple structures, to confirm adjustment on patient's implants and their positioning in the mouth. Provisional fixation with passive fit in the mouth and delivery back to the lab.



**8.**  
LAB\_ Final adjustments in the structure, where appropriate, aesthetic restoration, polishing and finishing. Reassembly of the replicas according to the actual relation in the mouth and cementing of both parts. Aesthetic restoration, polishing and finishing.



**9.**  
CLINIC\_ Placement of the prosthesis to the patient.

## 5. COMPATIBILITY CHART

| TRADEMARK            | PAGES | TRADEMARK           | PAGES |
|----------------------|-------|---------------------|-------|
| AB IMPLANTS          | 8     | IMPLANT DIRECT      | 22    |
| ACE                  | 8     | JDENTAL CARE        | 22    |
| ADIN                 | 8     | KEYSTONE DENTAL     | 23    |
| ALPHABIO             | 9     | KLOCKNER            | 23    |
| ALPHADENT            | 9     | LEADER              | 24    |
| ANCLADEN             | 9     | MEDENTIS            | 24    |
| ANTHOGYR             | 10    | MEGAGEN             | 24    |
| ARDS                 | 10    | MICRODENT           | 25    |
| ASTRA                | 10    | MIS                 | 26    |
| AVINENT              | 11    | MOZO GRAU           | 26    |
| BEGO                 | 11    | MPI                 | 27    |
| BIOGENESIS           | 12    | NEOBIOTECH          | 27    |
| BIOHORIZONS          | 13    | NEODENT             | 28    |
| BIOMET 3i            | 13    | NOBEL               | 28    |
| BIONER               | 13    | NT TRADING          | 29    |
| BTI                  | 14    | OSSTEM IMPLANT      | 29    |
| BTK                  | 14    | OSTEOPLUS           | 29    |
| B&W                  | 15    | OXTEIN              | 30    |
| CAMLOG               | 16    | PALTOP              | 30    |
| CON. SISTEMA PROTESE | 16    | PHIBO               | 31    |
| CORTEX               | 17    | RADHEX IMPLANTS     | 31    |
| COWELLMEDI           | 17    | SGS                 | 31    |
| DENTALTECH           | 17    | SIGNO VINCES        | 32    |
| DENTIS               | 18    | STRAUMANN           | 32    |
| DENTIUM              | 18    | SWEDEN & MARTINA    | 33    |
| DENTSPLY SIRONA      | 18    | TRI DENTAL IMPLANTS | 33    |
| DIO IMPLANTS         | 19    | TRINON              | 33    |
| EDISON MEDICAL       | 19    | ZIACOM (OSSEOLIFE)  | 34    |
| ELITE MEDICA         | 19    | ZIMMER              | 34    |
| EUROTEKNIKA          | 19    |                     |       |
| GALIMPLANT           | 20    |                     |       |
| GMI                  | 20    |                     |       |
| GT MEDICAL           | 21    |                     |       |
| HI TECH              | 21    |                     |       |



**[www.ballheadtechnology.com](http://www.ballheadtechnology.com)**

**(+34) 973 831 391 · [info@ballheadtechnology.com](mailto:info@ballheadtechnology.com)**